






# False Bottom Tubes






Wide range of tube options tailored for laboratory automation









- **NEW!** 13 x 75 mm
- Conical false bottom especially suitable for small volumes and platform compatibility
- Rounded base to accommodate rack requirements
- Fulfills IATA requirements for sample transport
- Suitable for deep freezing
- Centrifugation to 4,500 x g, 10 min.
- Convenient, compact StackPack packaging

# Ordering information • False bottom tubes ø 13 and 15 mm and compatible caps

Order No.	Length/Ø	Volume	Version	Tube base	Packaging
60.614.010	75/13 mm	2.5 ml	without cap, bulk packaging	round	1,000 pcs./box
60.614.015	75/13 mm	2.5 ml	without cap, StackPack packaging	round	2,000 pcs./box
<b>65.647</b>	<b>65.647.020</b>	<b>65.806</b>	<b>65.1379</b>	<b>65.163</b>	
Archiving cap for Ø 13 to 16 mm tubes	Archiving cap with filter for Ø 13 to 16 mm tubes	Push cap for Ø 13 mm tubes	Universal cap for Ø 12 to 16 mm tubes (additional colours available on request)	Screw cap Ø 13 mm for false bottom tubes	
					
5,000 pcs./box	5,000 pcs./box	10,000 pcs./box	10,000 pcs./box	5,000 pcs./box	

Order No.	Length/Ø	Volume	Version	Tube base	Packaging
60.617.010	92/13 mm	3.5 ml	without cap, without print	round	1,000 pcs./box
62.617	92/13 mm	3.5 ml	screw cap assembled, with print	round	
60.613.010	92/13 mm	3.5 ml	without cap, without print	flat	
62.613	92/13 mm	3.5 ml	screw cap assembled, with print	flat	
<b>65.647</b>	<b>65.647.020</b>	<b>65.806</b>	<b>65.1379</b>	<b>65.728</b>	
Archiving cap for Ø 13 to 16 mm tubes	Archiving cap with filter for Ø 13 to 16 mm tubes	Push cap for Ø 13 mm tubes	Universal cap for Ø 12 to 16 mm tubes (additional colours available on request)	Screw cap for S-Monovette® and Ø 13 mm false bottom tubes	
					
5,000 pcs./box	5,000 pcs./box	10,000 pcs./box	10,000 pcs./box	5,000 pcs./box	

Order No.	Length/Ø	Volume	Version	Tube base	Packaging
60.612.010	92/15.3 mm	5.0 ml	without cap, without print	round	1,000 pcs./box
62.612	92/15.3 mm	5.0 ml	screw cap assembled, with print	round	
60.611.011	92/15.3 mm	5.0 ml	without cap, without print, amber	flat	
60.611.010	92/15.3 mm	5.0 ml	without cap, without print	flat	
62.611	92/15.3 mm	5.0 ml	screw cap assembled, with print	flat	
<b>65.647</b>	<b>65.647.020</b>	<b>65.803</b>	<b>65.1379</b>	<b>65.645 *</b>	<b>65.729</b>
Archiving cap for Ø 13 to 16 mm tubes	Archiving cap with filter for Ø 13 to 16 mm tubes	Push cap for Ø 15.7 mm tubes	Universal cap for Ø 12 to 16 mm tubes (additional colours available on request)	Screw cap for S-Monovette® and Ø 15 mm false bottom tubes	Screw cap for S-Monovette® Ø 15 mm and 15.3 mm tubes
					
5,000 pcs./box	5,000 pcs./box	5,000 pcs./box	10,000 pcs./box	5,000 pcs./box	5,000 pcs./box

\* recommended for RC 1200 from Sarstedt

SARSTEDT AG & Co.  
P.O. Box 12 20  
D-51582 Nümbrecht  
Phone +49 2293 3050  
Fax +49 2293 305-122  
info@sarstedt.com  
www.sarstedt.com

 **SARSTEDT**