







False Bottom Tubes






Wide range of tube options tailored for laboratory automation









- Conical false bottom especially suitable for small volumes and platform compatibility
- Rounded and flat base options
- Fullfills IATA requirements for sample transport
- Suitable for deep freezing
- Centrifugation to 4,500 x g, 10 min. (13 x 75 mm tubes)
- Compact StackPack packaging available for 13 x 75 mm tubes

Ordering information • False bottom tubes ø 13 and 15 mm and compatible caps

Order No.	Length/Ø	Volume	Version	Tube base	Packaging
60.614.010	75/13 mm	2.5 ml	without cap, bulk packaging	round	1.000 pcs./box
60.614.015	75/13 mm	2.5 ml	without cap, StackPack packaging	round	2.000 pcs./box
65.647	65.647.020	65.806.300	65.1379	65.163	
Universal anti- evaporation cap for Ø 13 to 16 mm tubes	Universal anti- evaporation cap with filter for Ø 13 to 16 mm tubes	Push cap for Ø 13 mm tubes	Universal cap for Ø 12 to 16 mm tubes (additional colors available on request)	Screw cap Ø 13 mm for false bottom tubes	
					
5,000 pcs./box	5,000 pcs./box	10,000 pcs./box	10,000 pcs./box	5,000 pcs./box	

Order No.	Length/Ø	Volume	Version	Tube base	Packaging
60.617.010	92/13 mm	3.5 ml	without cap, without print	round	1.000 pcs./box
62.617	92/13 mm	3.5 ml	screw cap assembled, with print	round	
60.613.010	92/13 mm	3.5 ml	without cap, without print	flat	
62.613	92/13 mm	3.5 ml	screw cap assembled, with print	flat	
65.647	65.647.020	65.806.300	65.1379	65.728	
Universal anti- evaporation cap for Ø 13 to 16 mm tubes	Universal anti- evaporation cap with filter for Ø 13 to 16 mm tubes	Push cap for Ø 13 mm tubes	Universal cap for Ø 12 to 16 mm tubes (additional colors available on request)	Screw cap for S-Monovette® and Ø 13 mm false bottom tubes	
					
5,000 pcs./box	5,000 pcs./box	10,000 pcs./box	10,000 pcs./box	5,000 pcs./box	

Order No.	Length/Ø	Volume	Version	Tube base	Packaging
60.612.310	92/15.3 mm	5.0 ml	without cap, without print	round	1.000 pcs./box
62.612.300	92/15.3 mm	5.0 ml	screw cap assembled, with print	round	
60.611.311	92/15.3 mm	5.0 ml	without cap, without print, amber	flat	
60.611.310	92/15.3 mm	5.0 ml	without cap, without print	flat	
62.611.300	92/15.3 mm	5.0 ml	screw cap assembled, with print	flat	
65.647	65.647.020	65.803.300	65.1379	65.645 *	65.729.300
Universal anti- evaporation cap for Ø 13 to 16 mm tubes	Universal anti- evaporation cap with filter for Ø 13 to 16 mm tubes	Push cap for Ø 15.7 mm tubes	Universal cap for Ø 12 to 16 mm tubes (additional colors available on request)	Screw cap for S-Monovette® and Ø 15 mm false bottom tubes	Screw cap for S-Monovette® Ø 15 mm and 15.3 mm tubes
					
5,000 pcs./box	5,000 pcs./box	5,000 pcs./box	10,000 pcs./box	5,000 pcs./box	5,000 pcs./box

* recommended for RC 1200 from Sarstedt

SARSTEDT Inc.
1025 St. James Church Road
P.O. Box 468
Newton, NC 28658-0468
Tel. (800) 257-5101
Fax (828) 465-4003
customerservice.us@sarstedt.com
www.sarstedt.com

